

**BLASTOCYST FORM**

<b>Experiment no.</b>		<b>Date</b>		<b>Time</b>		<b>Initial</b>	
-----------------------	--	-------------	--	-------------	--	----------------	--

Day 8	Well	Date		Kinetic score			Morphology score		
		Total COC	TOTAL BL	BL. Rates %	1 BL	2 XB	3 H	1 Poor	2 Good
	1-1								
	1-2								
	1-3								
	1-4								
	2-1								
	2-2								
	2-3								
	2-4								
	3-1								
	3-2								
	3-3								
	3-4								
	4-1								
	4-2								
	4-3								
	4-4								

Calculation of scores: # BL x 1 + # XB x 2 + # H x 3 = A

Total: # BL + # XB + # H = B

Score: A/B

BL = Blastocyst

XB = Expanded Blastocyst

H = Hatched