

**BLASTOCYST FORM**

<b>Experiment no.</b>		<b>Date</b>		<b>Time</b>		<b>Initial</b>	
-----------------------	--	-------------	--	-------------	--	----------------	--

Day 8	Total COC	Date		Kinetic score			Morphology score		
		TOTAL BL	BL. Rates %	1 BL	2 XB	3 H	1 Poor	2 Good	3 Excellent
1-1									
1-2									
1-3									
1-4									
2-1									
2-2									
2-3									
2-4									
3-1									
3-2									
3-3									
3-4									
4-1									
4-2									
4-3									
4-4									

**Calculator of scores:**  
 Kinetics calculation: #BL x 1 + #XB x 2 + #H x 3 = K  
 Morphology calculation: #Poor x 1 + #Good x 2 + #Excellent x 3 = M  
 Total Blastocyst calculation: #BL + #XB + #H = T  
 Kinetics score: K/T  
 Morphology score: M/T

BL = Blastocyst  
 XB = Expanded Blastocyst  
 H = Hatched